


This drawing includes confidential information which is the property of ASHCROFT INC. This copy is loaned for a specific purpose, with the understanding that it will not be reproduced or used for any other purpose or disclosed to others, and it will be returned on demand.

ASHCROFT[®]

250 EAST MAIN STREET
STRATFORD, CT 06614



NOTES:

- 1) Process Fitting is G1/4 B Male
- 2) Instrument Fitting is G1/4 Female
- 3) Side closest to pressure tube loop is the "Instrument Side"

REV.	ECO	DESCRIPTION	BY	DATE
A	9802	INITIAL RELEASE	BEN WESLEY	3/1/2018

This print certified correct for		General Dimensions		
Customer's Name		G 1/4 Male-to-Female		
		Product Code: 13 2198 XX 27		
		Series 2198 Capillary Siphon		
Customer's Order No.		Drawn	Date	Date of Issue
Ashcroft Order No.		Ben Wesley	16-Jan-18	SHEET-1
Certified by	Date	070A4300		A